



ARNOLD PALMER
HOSPITAL
For Children

ECMO Initiation: Technical and Physiologic Considerations

Presented by
Jeffrey Amendola CCP CES-P

Chief Pediatric Perfusionist, Co-Medical Director APMC ECMO

February 20, 2025

Disclosures:

- I have no disclosures

Pre-Initiation Considerations

MDR

- Multi-disciplinary rounds/call
- Determine candidacy and staffing/equipment needs

VA or VV

Patient weight

- Circuit size
- Right size circuit for lowest prime and lowest surface area
- 1/4 circuit up to 24kg
- 3/8 circuit 25 kg and above
- Hybrid line setup
- Shunt for access
- We don't use a Bridge

Cannula selection

- Cannulation site
- Size appropriate for required flow
 - Sepsis
 - Low cardiac output

Checklist

Circuit Prime Considerations

Circuit Prime Considerations

Blood or Crystalloid

- Whole blood emergency access
 - Calcium and clotting considerations
 - HCT and volume concerns
- Washed blood
 - Time delay
- Unwashed RBC
 - Storage lesions
 - Potassium levels

Balanced prime

- CO₂ in neonates
Respiratory disease
- Congenital Diaphragmatic Hernia(CDH)
 - Neuro protective
- Oxygen level appropriate
- Terumo CDI inline blood gas analyzer

Circuit Prime Considerations

Other prime constituents

- Albumin
- Sodium Bicarbonate
- Calcium
- Heparin?

Hemofilter

- Wash blood prime
- Hemo-concentrate prime volume

Cannulation Considerations

Cannulation Considerations

Cannulas are selected and appropriate

- Extra supplies
- Cannula
 - Size up or down if anatomically challenging
- Connectors
- Tubing
- Scissors

Heparin

- Patient condition
- Bleeding, post surgical
- Low bleed risk
 - 50-100u/kg loading
 - 100-200u prime

Gas flow

- Blender settings appropriate for ECMO modality
- Plugged into wall
 - Not on tank
- Electric plugged in to RED
- Heater Cooler temp
- Shunt to recirculate prime during cannulation

Initiation of flow considerations

Initiation of flow considerations

Clamps off

- Establish target flow slowly
- 30-60 seconds

Initiation of flow considerations

Difficulties establishing flow: Circuit

- Obstructions to flow/drainage
- Alarms
 - High Venous inlet pressure
 - P1 P2 pressures, Delta P
- Clamps
 - Surgeon or you
- Bridges/Shunts open/closed affecting flow
 - Flow probes pre/post shunt
- RPM in a centrifugal pump
 - High enough to establish forward flow
- Roller pump
 - Servo-regulator
- Cannula
 - Placement/depth
 - Size
 - Anatomy
 - Patient positioning
- Tubing kinks

Initiation parameters: Patient VA

Initiation parameters: Patient VA

Blood pressure Adequate

- Flow requirements met?
 - Lactate load requires higher flow to clear Lactate
 - VA Target Flow of 100ml/kg may not be enough
 - Inotropes weaning or needed
 - Intra operative failure to wean from CPB

Arterial Sats

- Pre and Post Ductal sats
 - Higher flows required?

SVO2 appropriate range

Arterial line:Venous line color change?

Initiation parameters: Patient VA

Heart Rate

- Cardiac Stun
 - Ischemia
 - Calcium low
 - Potassium high
 - Need for a vent or decompression

Hemodilution

- Low HGB
 - DO₂ monitored

Fluid overload

- Over resuscitation
- May need to pull volume off circuit to regain pulsatility on Art Line

Initiation Parameters: VV

Initiation Parameters: VV

Recirculation

- Dual Lumen cannula position
 - Increasing flows can increase the problem, not fix them
 - Access and return cannulas in close proximity in traditional 2 cannula VV configurations
- Saturations
 - Arterial sats achieved?
May need to accept lower initial levels
 - Venous saturations unusually high indicate recirculation
- Heart Rate
 - Cardiac function
 - Essential for VV support
 - Arrhythmia
 - May need to convert to VA
 - Do you have cannulas to switch modes in an emergency
- Ventilation settings
 - Maximum vent settings adjusted to Safe vent settings

Special populations and true emergencies

Special populations and true emergencies

Palliative shunts

- Central shunt should be clipped to avoid overperfusion of the pulmonary system
- Coronary steal/ diastolic runoff to the pulmonary system
 - Compromising coronary perfusion

Cavo-pulmonary anastomosis(Glenn)

- Neck cannulation not advised due to difficult drainage/flow issues from anatomy
- Bi-Directional Bi-Lateral Glenn?

Fontan circulation

- Central cannulation in fresh procedure
- May need extra venous limb to drain VV-A
- Older failing Fontan vs a fresh extra-cardiac

Impella or other LVAD

- Impella can be used to offload LV if currently on support.
- IABP

Special population and True Emergencies

Special population and True Emergencies

Air introduction into circuit

Source

Patient

- Lines not secured or improper cannula connection at sterile field
- Accessing a venous port with high venous line negative pressure

Circuit

- Connection failure
- Accessing or Turning a stopcock incorrectly
 - Need for non occlusive caps on all port
- Infusions to venous side of circuit
- Avoid Arterial line

Air Embolism to patient

Cannula dislodgment/cannulation site failure

- Improper cannula suturing
- Frail tissue tear/failure
- Peripheral artery/vein tear or upon introduction

Special population and True Emergencies

Failure to Oxygenate

- Gas source
- O2 line disconnection
- Membrane Failure

Pump failure

- Decoupling/improper placement in motor housing
- Mechanical console failure
- Pump boot rupture


Bleeding

- Intra/post operative
 - No anticoagulation
 - MTP may be required to maintain circuit flow

Thank You

Establish a system you have
confidence in and rely on it when
you get into tough situations.

Arnold Palmer

 quote fancy



Great moments are born from Great Opportunity- Herb Brooks
February 22, 1980

References

Main section goes here and the copy box is as wide as the slide. Please do not increase the point size of any of the text. Main section goes here and the copy box is as wide as the slide.

- This is where the first copy point can go.
- Body copy on our presentation slides is ALWAYS 18 point.
- If bullet points are complete sentences, you should use punctuation on each one.
- You can also create a more graphic slide and put the details down into the notes area.
- This slide is nice to help break up your presentation and could just have photos instead of so much copy content.
- The first-indent level of bullet points is ALWAYS an open orange circle.



ECMO Initiation

Presented by
Andrew J. Snyder, Senior Vice President
Marketing & Communications

March 12, 2022

